

# Gen4 Sealed Call Cord



Why bother designing a new call cord? What could be improved on the simplest nurse signaling device? Everything. The design of the new Gen4 Sealed Call Cord was born from market demands, customer suggestions, and state-of-the-art design technology. It incorporates Curbell's 45 years of experience as the market leader to redefine the call cord.

## INFECTION CONTROL

- ▶ Proprietary sealed design keeps infectious contaminants out of the pendant cavity
- ▶ Smooth shape and minimal creases allow for easy, comprehensive cleaning, and reduce the potential for cross-contamination

## PATIENT SATISFACTION

- ▶ Modern, comfortable, ergonomic design fits naturally in the hand
- ▶ Large button is easy for patients to locate and operate
- ▶ Call assurance LEDs for supporting systems

## OXYGEN SAFETY

- ▶ UL Listing involves testing with the nurse call system and verifying for safe use in an oxygen-enriched atmosphere
- ▶ Visit [www.oxygensafe.com](http://www.oxygensafe.com) for more information

## REGULATORY COMPLIANCE

According to the Food and Drug Administration, nurse call systems and their accessories are Class II Medical Devices (21 CFR 890.3710). As a manufacturer of these devices, Curbell is required to meet specific requirements, including:

- ▶ Compliance with the current good manufacturing practices (cGMP) identified in the FDA's Quality System Regulation (21 CFR 820), and the Medical Device Reporting Regulation
- ▶ Registration of our manufacturing facility with the FDA
- ▶ Listing of the device with the FDA